


目录 | CATALOG
ENTERPRISE DIRECTORY

- ▶ 医疗床系列
Hospital Beds Series
- ▶ 居家/养老院护理系列
Home/Nursing Home Nursing Series
- ▶ 轮椅系列
Wheelchair Series
- ▶ 急救系列
First aid Series
- ▶ 配件
Accessories
- ▶ 配套设施
Supporting Facilities

ELECTRIC 医疗床系列 APITAL
Medical bed series



Electric Hospital Bed

CH-C03

外部长度 External Length	2150mm
外部宽度 External Width	1030mm
床面高度 Bed Height	550-780mm
床垫尺寸 Mattress Size	1960*900mm
安全承重 Safe Load	150kg

遥控器 Remote Controller	1pc	引流钩 Drainage Hook	2pcs
中控脚轮 Central Control Casters	4pcs	输液杆 IV Pole	1set
内嵌控制按键的护栏 Side rail with integrated control panel	1set	输液杆插孔 IV Pole Prevision	2sets
床头/床尾 Head/Footboard	1set		
电机 Linear Motor	4pcs		

可选配置 Optional

康复吊环 Monkey Pole	<input type="checkbox"/>
防褥疮气垫 Air Mattress	<input type="checkbox"/>
床垫 Mattress	<input type="checkbox"/>
杂物架 Shelf	<input type="checkbox"/>
ABS床面 ABS Bed Surface	<input type="checkbox"/>
后背电池 Backup Battery	<input type="checkbox"/>
机械CPR Manual CPR	<input type="checkbox"/>
电子CPR Electrical CPR	<input type="checkbox"/>
防撞轮 Wall Bumper Wheel	<input type="checkbox"/>

品名: Electric hospital bed

电动医疗床 CH-C01

产品参数 (Product parameter)

床面长度 External Length	1960±20mm
床面宽度 External Width	900±20mm
床面高度 Bed Height	420/700mm
背部倾斜 Back Tilt	≥ 75° ± 10°
腿部倾斜 Leg tilt	≥ 35° ± 10°

产品介绍 Product Description

- 1、床体采用冷轧钢板冲压而成；
- 2、床头床尾板采用高端ABS注塑材质，外形美观，装卸自如；
- 3、四片高品质缓降护栏，铝合金支架，方便操作，结实耐用；
- 4、高品质万向静音脚轮，中控制车、移动灵活、稳定性好；
- 5、超静音直流线性推杆。

1. High quality engineering ABS plastic head & foot board.
2. Overall punching cold carbon steel bed board which is stable and reliable.
3. ABS plastic hydraulic side rails to protect patient from falling off bed.
4. ABS double sides central locking control mute casters help to move the bed easily.
5. ROHS & CE certificated high quality and powerful electric pushrods and controllers.

功能介绍 Function Description



Electric hospital bed 电动医疗床

CH-C02 尺寸 (Size) 2150*1030*420/700mm



Electric hospital bed 电动医疗床

CH-E05A 尺寸 (Size) 2150*1030*420/700mm



Electric hospital bed 电动医疗床

CH-C04 尺寸 (Size) 2150*1030*420/700mm



Electric hospital bed 电动医疗床

CH-E05B 尺寸 (Size) 2150*1030*420/700mm



Electric hospital bed 电动医疗床

CH-E05C 尺寸 (Size) 2150*960*420/700mm



Electric hospital bed 电动医疗床

CH-E03A 尺寸 (Size) 2150*1030*420/700mm



Electric hospital bed 电动医疗床

CH-E03C 尺寸 (Size) 2150*960*420/700mm




Electric hospital bed 电动医疗床

CH-E03E 尺寸 (Size) 2150*1030*420/700mm



Electric hospital bed 电动医疗床

CH-E05D 尺寸 (Size) 2150*960*420/700mm



Electric hospital bed 电动医疗床

CH-E03B 尺寸 (Size) 2150*1030*420/700mm



Electric hospital bed 电动医疗床


CH-E03D 尺寸 (Size) 2150*960*420/700mm



Electric hospital bed 电动医疗床

CH-E03F 尺寸 (Size) 2150*1030*420/700mm





Manual Hospital Bed

CH-M04A

外部长度 External Length	2150mm
外部宽度 External Width	1030mm
床面高度 Bed Height	420-700mm
床垫尺寸 Mattress Size	1960*900mm
安全承重 Safe Load	250kg

摇柄 Manual Crank	4pc	输液杆 IV Pole	1set
脚轮 5 inch Casters	4pcs	输液杆插孔 IV Pole Prevision	2sets
护栏 Foldable Side Rails	1set		
床头/床尾 Head/Footboard	1set		
引流沟 Drainage Hook	2pcs		

可选配置 Optional

康复吊环 Monkey Pole	<input type="checkbox"/>
防褥疮气垫 Air Mattress	<input type="checkbox"/>
床垫 Mattress	<input type="checkbox"/>
杂物架 Shelf	<input type="checkbox"/>
中控脚轮 Central Control Casters	<input type="checkbox"/>
防撞轮 Wall Bumper Wheel	<input type="checkbox"/>

品名: Manual Orthopedic bed

手动医疗床 CH-M06

产品参数 (Product parameter)

床面长度 External Length	1960±20mm
床面宽度 External Width	900±20mm
床面高度 Bed Height	420/700mm
背部倾斜 Back Tilt	≥ 75° ± 10°
腿部倾斜 Leg tilt	≥ 35° ± 10°

产品介绍 Product Description

1. 高品质ABS复合床头板带防撞轮。
 2. 冷碳钢管床架，腿部分体设计。
 3. 高品质的铝合金6档可折叠护栏。
 4. 带独立锁定装置的ABS静音脚轮。
 5. 高品质滚珠丝杠和可折叠不锈钢摇柄。
1. Head board and feet board are made of high quality ABS engineering plastic.
 2. Cold carbon steel welded bed frame which is stable and reliable.
 3. High-quality aluminum alloy foldable side rail.
 4. ABS silent caster with independent locking device.
 5. Equipped with 60*60mm high quality pushing screws with stainless steel foldable cranks.

功能介绍 Function Description



Manual hospital bed 手摇床

CH-M03A 尺寸 (Size) 2150*960*420/700mm




Manual hospital bed 手动医疗床

CH-M03B 尺寸 (Size) 2070*960*420/700mm




Manual hospital bed 手摇床

CH-M03C 尺寸 (Size) 2150*960*420/700mm



Manual hospital bed 手摇床

CH-M03D 尺寸 (Size) 2150*1030*420/700mm



品名: Manual hospital bed

手动医疗床 CH-M02E

产品参数 (Product parameter)

床面长度 External Length	1960±20mm
床面宽度 External Width	900±20mm
床面高度 Bed Height	500±20mm
背部倾斜 Back Tilt	≥ 75° ± 10°
腿部倾斜 Leg tilt	≥ 35° ± 10°

产品介绍 Product Description

1. 高品质ABS复合床头板;
2. 冷碳钢床架, 整体冲孔床面;
3. PP缓冲护栏, 铝合金支架;
4. 带中控制车的高品质静音脚轮;
5. 高品质滚珠丝杠和ABS可折叠摇柄。


1. Head board and feet board are made of high quality ABS engineering plastic.
2. Cold carbon steel welded bed frame which is stable and reliable.
3. ABS plastic hydraulic side rails to protect patient from falling off bed.
4. ABS mute casters help to move the bed easily.
5. Equipped with 60*60mm high quality pushing screws with stainless steel foldable cranks.

功能介绍 Function Description



Manual hospital bed 手摇床

CH-M02 尺寸 (Size) 2150*960*500mm




Manual hospital bed 手摇床

CH-M02A 尺寸 (Size) 2150*960*500mm



Manual hospital bed 手摇床

CH-M02C 尺寸 (Size) 2070*960*500mm



Manual hospital bed 手摇床

CH-M02D 尺寸 (Size) 2150*960*500mm



Manual hospital bed 手摇床

CH-M02B 尺寸 (Size) 2150*960*500mm



Manual hospital bed 手摇床

CH-M02F 尺寸 (Size) 2150*1030*500mm



Manual hospital bed 手摇床

CH-M01B 尺寸 (Size) 2070*960*500mm




Manual hospital bed 手摇床

CH-M01A 尺寸 (Size) 2070*960*500mm



ELECTRONIC NURSING
居家/养老院护理系列
Home/Nursing Home Nursing Series



Electric Nursing Bed

CN-E20

外部长度 External Length	2050mm
外部宽度 External Width	950mm
床面高度 Bed Height	500mm
床垫尺寸 Mattress Size	1960*900mm
安全承重 Safe Load	150kg

遥控器 Remote Controller	1pc	引流沟 Drainage Hook	2pcs
脚轮 5 inch Casters	4pcs	输液杆 IV Pole	1set
护栏 Foldable Side Rails	1set	输液杆插孔 IV Pole Prevision	2sets
床头/床尾 Head/Footboard	1set		
电机 Linear Motor	4pcs		

可选配置 Optional

康复吊环 Monkey Pole	<input type="checkbox"/>
防褥疮气垫 Air Mattress	<input type="checkbox"/>

品名: Electric Nursing bed

电动护理床 CN-E16

产品参数 (Product parameter)

床面长度 External Length	1960±20mm
床面宽度 External Width	900±20mm
床面高度 Bed Height	480±20mm
背部倾斜 Back Tilt	≥ 75° ± 10°
腿部抬起 Leg tilt	30° ± 5°
腿部降落 Leg landing	65° ± 10°
左翻/右翻 Turn Left/Right	55° ± 5°

产品介绍 Product Description

1. 木质床头板;
 2. 冷碳钢床架, 可折叠护栏;
 3. 带独立锁定装置的ABS双面静音脚轮;
 4. 高品质电动推杆和控制器;
 5. 电动便孔。
1. Wooden head & feet boards which apply to home use.
 2. Cold carbon steel welded bed frame which is stable and reliable.
 3. ABS silent casters help to move the bed easily.
 4. ROHS & CE certificated high quality and powerful electric pushrods and controllers.
 5. Electric push toilet with waterproof mattress which is convenient for patient to use and keep clean.

功能介绍 Function Description



Electric nursing bed 电动护理床

CN-E02 尺寸 (Size) 2070*960*480mm



Electric nursing bed 电动护理床

CN-E03 尺寸 (Size) 2070*960*480mm



Electric nursing bed 电动护理床

CN-E04 尺寸 (Size) 2070*1000*480mm



Electric nursing bed 电动护理床

CN-E15 尺寸 (Size) 2070*960*480mm



Electric nursing bed 电动护理床

CN-E01 尺寸 (Size) 2070*960*480mm



Electric nursing bed 电动护理床

CN-E12 尺寸 (Size) 2070*960*480mm




Electric nursing bed 电动护理床

CN-E06 尺寸 (Size) 2070*960*480mm



Electric back lifter 电动起背器

CN-E08 尺寸 (Size) 1900*900mm



品名: Manual Nursing bed

手动护理床 CN-M16

产品参数 (Product parameter)

床面长度 External Length	1960±20mm
床面宽度 External Width	900±20mm
床面高度 Bed Height	480±20mm
背部倾斜 Back Tilt	≥ 75° ± 10°
腿部抬起 Leg tilt	30° ± 5°
腿部降落 Leg landing	65° ± 10°
左翻/右翻 Turn Left/Right	55° ± 5°

产品介绍 Product Description

1. 高品质ABS床头板;
 2. 冷碳钢床架, PP缓冲护栏, 铝合金支架;
 3. 带独立锁定装置的高品质ABS静音脚轮;
 4. 高品质滚珠丝杠和不锈钢可折叠摇柄;
 5. 速开式便孔。
1. Head board and feet board are made of high quality ABS engineering plastic.
 2. Cold carbon steel welded bed frame which is stable and reliable.
 3. ABS mute casters help to move the bed easily.
 4. Equipped with 60*60mm high quality pushing screws with stainless steel foldable cranks.
 5. Crank style toilet can be quickly opened and closed.

功能介绍 Function Description



Manual Nursing Bed 手动护理床

CN-M15 尺寸 (Size) 2070*960*480mm



Manual Nursing Bed 手动护理床

CN-M10 尺寸 (Size) 2070*960*480mm



Manual Nursing Bed 手动护理床

CN-M12 尺寸 (Size) 2070*960*480mm



Manual Nursing Bed 手动护理床

CN-M15B 尺寸 (Size) 2070*960*480mm



Electric hospital bed 电动床

CN-E21 尺寸 (Size) 2150*1030*500mm



Manual hospital bed 手摇床

CN-M02B 尺寸 (Size) 2070*960*500mm



Manual hospital bed 手摇床

CN-M21 尺寸 (Size) 2150*1030*500mm




Manual hospital bed 手摇床

CN-M02A 尺寸 (Size) 2070*960*500mm



Electric hospital Bed 电动床

CN-E31 尺寸 (Size) 2070*960*420/700mm




Manual Hospital Bed 手摇床

CN-M31 尺寸 (Size) 2070*960*420/700mm



Manual Hospital Bed 手摇床

CN-M02C 尺寸 (Size) 2070*960*500mm



Manual Hospital Bed 手摇床

CN-M02D 尺寸 (Size) 2070*960*500mm



ELECTRI 轮椅系列 ELCHAIR
Wheelchair series

Electric Wheelchair

电动轮椅

PW-6019

整体尺寸 Overall size	118x63x93cm
框架材质 Frame	铝合金
前轮尺寸 Front Wheel	8寸
后轮尺寸 Rear Wheel	12寸
座宽 Seat Width	45cm
承重 Load Capacity	100kg
电池 Battery	12AH 24V
续航 Endurance	15-20km



功能介绍 >>

Width
座宽



45cm

Depth
座深



45cm

Speed
行驶速度



km/h
0-6km/h

Voyage
续航



15-20km

Front wheel
前轮



8"

Rear wheel
后轮



12"

Motor
电机



250Wx2

Battery
电池



12AH 24V

PW-6016

电动轮椅车
Electric Wheelchair

整体尺寸 Overall size	118x63x93cm
框架材质 Frame	铝合金
前轮尺寸 Front Wheel	8寸
后轮尺寸 Rear Wheel	12寸
座宽 Seat Width	45cm
承重 Load Capacity	100kg
电池 Battery	12AH 24V
续航 Endurance	15-20km



PW-6001A

电动轮椅车
Electric Wheelchair

整体尺寸 Overall size	118x63x93cm
框架材质 Frame	碳钢
前轮尺寸 Front Wheel	10寸
后轮尺寸 Rear Wheel	16寸
座宽 Seat Width	45cm
承重 Load Capacity	100kg
电池 Battery	12AH 24V
续航 Endurance	15-20km



PW-6005A

电动轮椅车
Electric Wheelchair

整体尺寸 Overall size	118x63x93cm
框架材质 Frame	碳钢
前轮尺寸 Front Wheel	10寸
后轮尺寸 Rear Wheel	16寸
座宽 Seat Width	45cm
承重 Load Capacity	100kg
电池 Battery	12AH 24V
续航 Endurance	15-20km



PW-6002

电动轮椅车
Electric Wheelchair

整体尺寸 Overall size	118x63x93cm
框架材质 Frame	碳钢
前轮尺寸 Front Wheel	10寸
后轮尺寸 Rear Wheel	22寸
座宽 Seat Width	45cm
承重 Load Capacity	100kg
电池 Battery	12AH 24V
续航 Endurance	15-20km



Manual Wheelchair

手动轮椅

PW-8030

整体尺寸 Overall size	105x68x121cm
框架材质 Frame	碳 钢
前轮尺寸 Front Wheel	8寸
后轮尺寸 Rear Wheel	12寸
座 宽 Seat Width	45cm
坐 深 Seat Deep	41cm
承 重 Load Deep	100kg
包装尺寸 Package Size	93x26x97cm



功能介绍 >>



PW-8001

手动轮椅车
Manual Wheelchair

整体尺寸 Overall size	103x66x89cm
框架材质 Frame	碳 钢
前轮尺寸 Front Wheel	7寸
后轮尺寸 Rear Wheel	24寸
座 宽 Seat Width	45cm
座 深 Seat Deep	41cm
承 重 Load Capacity	100kg
包装尺寸 Package Size	93x26x97cm



PW-8003

手动轮椅车
Manual Wheelchair

整体尺寸 Overall size	103x66x89cm
框架材质 Frame	碳 钢
前轮尺寸 Front Wheel	7寸
后轮尺寸 Rear Wheel	24寸
座 宽 Seat Width	45cm
座 深 Seat Deep	41cm
承 重 Load Capacity	100kg
包装尺寸 Package Size	93x26x97cm



PW-8002

手动轮椅车
Manual Wheelchair

整体尺寸 Overall size	102x68x91cm
框架材质 Frame	碳 钢
前轮尺寸 Front Wheel	7寸
后轮尺寸 Rear Wheel	24寸
座 宽 Seat Width	45cm
座 深 Seat Deep	41cm
承 重 Load Capacity	100kg
包装尺寸 Package Size	93x26x97cm



PW-8004

手动轮椅车
Manual Wheelchair

整体尺寸 Overall size	103x66x89cm
框架材质 Frame	碳 钢
前轮尺寸 Front Wheel	7寸
后轮尺寸 Rear Wheel	24寸
座 宽 Seat Width	45cm
座 深 Seat Deep	41cm
承 重 Load Capacity	100kg
包装尺寸 Package Size	93x26x97cm



PW-8011

手动轮椅车
Manual Wheelchair

整体尺寸 Overall size	82*61*91cm
框架材质 Frame	碳钢
前轮尺寸 Front Wheel	6寸
后轮尺寸 Rear Wheel	16寸
座宽 Seat Width	41cm
座深 Seat Deep	41cm
承重 Load Capacity	100kg
包装尺寸 Package Size	75*28*70cm

**PW-8012**

手动轮椅车
Manual Wheelchair

整体尺寸 Overall size	85*61*91cm
框架材质 Frame	碳钢
前轮尺寸 Front Wheel	6寸
后轮尺寸 Rear Wheel	20寸
座宽 Seat Width	41cm
座深 Seat Deep	41cm
承重 Load Capacity	100kg
包装尺寸 Package Size	75*28*70cm

**PW-8020**

手动加宽轮椅车
Manual width Wheelchair

整体尺寸 Overall size	106*73*97cm
框架材质 Frame	碳钢喷涂
前轮尺寸 Front Wheel	7寸
后轮尺寸 Rear Wheel	24寸
座宽 Seat Width	49cm
座深 Seat Deep	45cm
承重 Load Capacity	100kg
包装尺寸 Package Size	93*26*97cm

**PW-8007**

手动轮椅车
Manual Wheelchair

整体尺寸 Overall size	102*68*91cm
框架材质 Frame	碳钢
前轮尺寸 Front Wheel	7寸
后轮尺寸 Rear Wheel	24寸
座宽 Seat Width	45cm
座深 Seat Deep	41cm
承重 Load Capacity	100kg
包装尺寸 Package Size	93*26*97cm

**PW-8005**

手动轮椅车
Manual Wheelchair

整体尺寸 Overall size	102*68*91cm
框架材质 Frame	碳钢
前轮尺寸 Front Wheel	7寸
后轮尺寸 Rear Wheel	24寸
座宽 Seat Width	45cm
座深 Seat Deep	41cm
承重 Load Capacity	100kg
包装尺寸 Package Size	93*26*97cm

**PW-8006**

手动轮椅车
Manual Wheelchair

整体尺寸 Overall size	102*68*91cm
框架材质 Frame	碳钢
前轮尺寸 Front Wheel	7寸
后轮尺寸 Rear Wheel	24寸
座宽 Seat Width	45cm
座深 Seat Deep	41cm
承重 Load Capacity	100kg
包装尺寸 Package Size	93*26*97cm

**PW-8008**

手动轮椅车
Manual Wheelchair

整体尺寸 Overall size	102*68*91cm
框架材质 Frame	碳钢喷涂
前轮尺寸 Front Wheel	7寸
后轮尺寸 Rear Wheel	24寸
座宽 Seat Width	45cm
座深 Seat Deep	41cm
承重 Load Capacity	100kg
包装尺寸 Package Size	93*26*97cm

**PW-8009**

手动轮椅车
Manual Wheelchair

整体尺寸 Overall size	102*68*91cm
框架材质 Frame	碳钢
前轮尺寸 Front Wheel	7寸
后轮尺寸 Rear Wheel	24寸
座宽 Seat Width	45cm
座深 Seat Deep	41cm
承重 Load Capacity	100kg
包装尺寸 Package Size	93*26*97cm




FIRST AID 急救系列 SERIES
First aid series

luxury transporter

PZ-OCT105


产品描述 (Product Description)

1. 框架采用优质碳钢焊接成型，对接轨道采用高强度铝合金型材；
 2. 车面及护栏采用进口 PP 工程塑料一次注塑成型；
 3. 两单体车可升降调节高度，两单车均采用高档中控刹车系统。
1. The frame is welded with carbon steel, and the butt rail is made of high-strength aluminum alloy;
 2. The bed surface and side rails are injection molded with high-quality PP engineering plastics;
 3. Both trolleys height can be adjusted, and casters are used central control brake systems.




主车尺寸	2080x720x780mm	摇柄	1pcs	引流沟	2pcs
副车尺寸	1960x740x660mm	方向轮	2pcs		
最大高度	610mm	带中控锁脚轮	2set		
最小高度	930mm	护栏	1set		


PZ-101 豪华转运车
Luxury Hydraulic Stretcher
尺寸 (Size) : 1980x800x610/930mm




PZ-102 豪华转运车
Luxury Hydraulic Stretcher
尺寸 (Size) : 1930x640x630/960mm



PZ-101 II 豪华转运车
Luxury Hydraulic Stretcher
尺寸 (Size) : 1885x800x630/960mm



PZ-OCT103 豪华转运车
Luxury Hydraulic Stretcher
尺寸 (Size) : 1930x640x500/900mm



Ambulance stretcher

救护车担架

CB-A1

产品描述 (Product Description)



- 1、床面采用泡沫垫制成，可调节背部使伤员感觉舒适。
- 2、机器非常方便：一个人足以将病人送到救护车担架。
- 3、机械折叠腿由担架两端的手柄控制。
- 4、当装入车辆时，腿可以自动折叠在担架下面。
- 5、在救护车上时，担架也可以通过固定装置锁定。
- 6、这种类型的担架的特点是轻便，耐用，易于消毒杀菌。

1. The bed surface is made of foam pad, which can adjust the back to make the wounded feel comfortable.
2. Machines are very convenient; One person is enough to carry the patient to the ambulance stretcher.
3. Mechanically folded legs are controlled by handles at each end of the stretcher.
4. Legs can automatically fold under the stretcher when loaded into the vehicle.
5. The stretcher can also be locked with a fixture while in an ambulance.
6. This type of stretcher is characterized as lightweight, durable, and easy to sterilize.



最高位 Top position	190*55*92cm
最低位 Lowest position	193*55*26cm
背部升降 Back Up	75° ± 5°
最大承重 Load Capacity	≤159kg

功能介绍



CB-A4 救护车担架

Ambulance stretcher

尺寸 (Size) : 2000x560x1000mm



CB-A2-2 救护车担架

Ambulance stretcher

尺寸 (Size) : 1910x570x910mm



CB-A3 救护车担架

Ambulance stretcher



尺寸 (Size) : 1950x570x380/910/870/780mm



CB-A1-4 椅式担架

Chair stretcher

尺寸 (Size) : 1960x550x260/980mm



CB-L1 楼梯担架
Stair stretcher
尺寸 (Size) : 900x510x910mm



CB-L1-1 楼梯担架
Stair stretcher
尺寸 (Size) : 960x550x260/980mm



CB-D3-3 脊椎板
Vertebral
尺寸 (Size) : 1830x450x50mm



CB-D3-4 脊椎板
Vertebral
尺寸 (Size) : 1870x500x40mm



CB-C1 铲式担架
Shovel stretcher
尺寸 (Size) : 2100x440x600mm



CB-C1-1 铲式担架
Shovel stretcher
尺寸 (Size) : 2140x420x70mm



CB-C2-1 铲式担架
Shovel stretcher
尺寸 (Size) : 1723x436x73mm



CB-D1-1 折叠担架
Fold stretcher
尺寸 (Size) : 2140x550x120mm



CB-D1-2 折叠担架
Fold stretcher
尺寸 (Size) : 2200x220x150mm



CB-C2 铲式担架
Shovel stretcher
尺寸 (Size) : 1620x450x70mm



CB-D2-3 折叠担架
Fold stretcher
尺寸 (Size) : 1850x500x240mm



CB-D2 折叠担架
Fold stretcher
尺寸 (Size) : 1850x500x30mm



CB-D2 折叠担架
Fold stretcher
尺寸 (Size) : 1850x500x20mm





ACCESSORIES

配件系列
Accessories Series

PS- AM01 充气床垫
Air mattress

尺寸/Size:1950x900x100mm



PS- AM02 充气床垫
Air mattress

尺寸/Size:1950x900x100mm



PS- AM04 充气床垫
Air mattress

尺寸/Size:1950x900x100mm



PS- AM05 充气床垫
Air mattress

尺寸/Size:1950x900x100mm



PS- AM06 充气床垫
Air mattress

尺寸/Size:1950x900x100mm



PS - O1 移动餐桌 Overbed Table
尺寸/Size: 900x400x740-1130mm



PS - O2 ABS移动餐桌 ABS Overbed Table
尺寸/Size: 820x450x840-1070mm



PS - O3 移动餐桌 Overbed Table
尺寸/Size: 800x400x740-1070mm



PS - O05A 餐桌板
Dining Table
尺寸/Size:950x400mm



PS - O05B 木质餐桌板
Wooden Dining Table
尺寸/Size:950x360mm



PS - O6 可伸缩餐桌板
Retractable Dining Table
尺寸/Size: (830-1010)x360mm



PS- G01A ABS床头柜
ABS Bedside Cabinet
尺寸/Size:480x480x760mm



PS- G01C ABS床头柜
ABS Bedside Cabinet
尺寸/Size:480x480x780mm



PS- G01F ABS床头柜
ABS Bedside Cabinet
尺寸/Size:480x480x880mm



PS- G02 ABS床头柜
ABS Bedside Cabinet
尺寸/Size:450x420x780mm



PS- G03 ABS床头柜
ABS Bedside Cabinet
尺寸/Size:490x480x780mm



PS- G01A ABS床头柜
ABS Bedside Cabinet
尺寸/Size:480x480x760mm



PS- G01A 钢塑床头柜
Steel-plastic Bedside Cabinet
尺寸/Size:510x510x710mm





PS- G05A 钢塑床头柜
Steel-plastic Bedside Cabinet
尺寸/Size:510x510x710mm



PS- G06 钢塑床头柜
Steel-plastic Bedside Cabinet
尺寸/Size:465x455x810mm





SUPPORTING FACILITIES

配套设施
Supporting facilities



PS- GAT03 豪华麻醉车
Luxury Anesthesia Trolley
尺寸/Size:615x560x1470mm



PS- AT101 麻醉车
Anesthesia Trolley
尺寸/Size:690x500x1550mm



PS- AT02 麻醉车
Anesthesia Trolley
尺寸/Size:1065x635x1600mm



PS- ET16 急救车
Emergency Trolley
尺寸/Size:690x500x990mm



PS- ET20 急救车
Emergency Trolley
尺寸/Size:690x500x990mm



PS- MT02 送药车
Medicine Trolley
尺寸/Size:630x530x900mm



PS- ET05 急救车
Emergency Trolley
尺寸/Size:680x450x940mm



PS- ET06 急救车
Emergency Trolley
尺寸/Size:630x440x850mm



PS- ET07 急救车
Emergency Trolley
尺寸/Size:605x455x990mm



PS- CT102 治疗车
Clinical Trolley
尺寸/Size:620x410x940mm



PS- CT103 治疗车
Clinical Trolley
尺寸/Size:710x460x840mm



PS- CT108 治疗车
Clinical Trolley
尺寸/Size:830x510x900mm



PS- ET08 急救车
Emergency Trolley
尺寸/Size:620x470x990mm



PS- ET09 急救车
Emergency Trolley
尺寸/Size:620x470x990mm



PS- ET09H 急救车
Emergency Trolley
尺寸/Size:620x470x990mm



PS- CT01 治疗车
Clinical Trolley
尺寸/Size:775x635x945mm



PS- CT08 治疗车
Clinical Trolley
尺寸/Size:690x450x980mm



PS- CT09 治疗车
Clinical Trolley
尺寸/Size:710x620x1010mm



PS- SST15 手杖凳
Cane With Stool



PS- SST17 腋下拐杖
Under Arm Crutches



PS- SST18-I 肘拐杖
Elbow Cane



PS- SST19-III 移动输液架
Infusion Stand



PS- SST10 助行器
Walker



PS- SST10-IV 助行器
Walker



PS- SST11-II 坐便椅
Stubbom Chair



PS- SST11-VI 坐便椅
Stubbom Chair



PS- SST15-I 手杖凳
Cane With Stool



PS- SST17-I 腋下拐杖
Under Arm Crutches



PS- SST19 移动输液架
Infusion Stand



PS- SST19-IV 移动输液架
Infusion Stand



PS- SST10-I 助行器
Walker



PS- SST10-V 助行器
Walker



PS- SST11-III 坐便椅
Stubbom Chair



PS- SST12 沐浴凳
Bath Chair



PS- SST16 拐杖
Cane



PS- SST17-II 腋下拐杖
Under Arm Crutches



PS- SST19-I 移动输液架
Infusion Stand



PS- SST19-V 移动输液架
Infusion Stand



PS- SST10-II 助行器
Walker



PS- SST11 坐便椅
Stubbom Chair



PS- SST11-IV 坐便椅
Stubbom Chair



PS- SST12-I 沐浴凳
Bath Chair



PS- SST16-I 拐杖
Cane



PS- SST18 肘拐杖
Elbow Cane



PS- SST19-II 移动输液架
Infusion Stand



PS- SST20 烤灯
Baked lamp



PS- SST10-III 助行器
Walker



PS- SST11-I 坐便椅
Stubbom Chair



PS- SST11-V 坐便椅
Stubbom Chair



PS- SST12-II 沐浴凳
Bath Chair



PS- SST01 陪护椅
Accompany Chair

尺寸/Size:720x460x900mm



PS- SST04 候诊椅
Waiting Chair



PS- SST07 天轨
Roof Track



PS- SST02 检查椅
Examination Chair

尺寸/Size:1680x600x800mm



PS- SST05 护士凳
Nurse Chair



PS- SST08 屏风
Screen



PS- SST03 输液椅
IV Chair

尺寸/Size:700x1100x1150mm



PS- SST06 检查床
Examination Bed



PS- SST09 设备带
Equipment Trunking