

BKAMED Press Plus & Press2 Plus

Pressure-Activated Disposable Safety Lancets

New double-spring design

Adequate blood volume, minimal pain



Color	Press Plus Reorde#	Color	Press2 Plus Reorde#	Gauge	Penetration Depth
	05-233015		05-243015	30G	1.5mm
	05-232818		05-242818	28G	1.8mm
	05-232618		05-242618	26G	1.8mm
	05-232122		05-242122	21G	2.2mm
	05-232118		05-242118	21G	1.8mm
	05-231720B		05-241720B	17GBlade	2.0mm

BKAMED ImPress Pro & ImPress



Pressure-Activated Disposable Safety Lancets

Recommended for blood sugar monitoring devices



Tactile Contact

Stable and accurate positioning to collect an adequate volume of blood



Color	Impress Pro Reorde#	Color	Impress Reorde#	Gauge	Penetration Depth
	05-193015		05-173015	30G	1.5mm
	05-192818		05-172818	28G	1.8mm
	05-192618		05-172618	26G	1.8mm
	05-192118		05-172118	21G	1.8mm
	05-192122		05-172122	21G	2.2mm

Side button-activated Disposable Safety Lancet



◇ Accurate positioning
Raised target

◇ Large anti-slip button
Shorter activation travel

◇ Optimized launch structure
Minimized pain for patients, secure handling for staff

Fixes the target position before firing

Relieves patient tension by minimizing pressure on the fingertip.



Color	Reorder#	Gauge	Penetration Depth
	05-113018	30G	1.8mm
	05-112818	28G	1.8mm
	05-112618	26G	1.8mm
	05-112118	21G	1.8mm
	05-111720B	17G Blade	2.0mm

Depth-adjustable, Pressure-activated lancing Disposable Safety Lancet



◇ Increased operational efficiency
One-twist uncapping

◇ Easy identification and error prevention
Usage status indicator

◇ Accurate positioning
Extended target port


Low / Medium / High blood sample volume

◇ Puncture depth 1.3mm / 1.8mm / 2.3mm
The most suitable penetration depth can be configured on the spot, according to the patient's skin thickness and capillary circulation, increasing the sampling success rate.




Color	Reorder#	Gauge	Penetration Depth
	05-1528A	28G	1.3/1.8/2.3mm
	05-1526A	26G	1.3/1.8/2.3mm
	05-1523A	23G	1.3/1.8/2.3mm

Rear button activated Disposable Safety Lancet Gamma-sterilized, assuring sterility for the patient



Safer T-shape design

Best fit for nurses needing gentle skin contact
Large anti-slip rear button providing a stable and comfortable pressing surface.



◇ Quick and straight puncture
Rear firing spring

◇ Reduced jitter during puncture
Front damping spring



◇ Usage status indicator



◇ Stable positioning


Comfortable double-spring action



Color	Reorder#	Gauge	Penetration Depth
	05-213015	30G	1.5mm
	05-212818	28G	1.8mm
	05-212618	26G	1.8mm
	05-212118	21G	1.8mm
	05-211720B	17G Blade	2.0mm


Depth-adjustable, Rear button activated Disposable Safety Lancet

Single SKU helps simplify inventory management



One Lancet offers three depth settings (1.3, 1.8, and 2.3mm) for various patient conditions

Simply select the depth that provides the greatest comfort—while at the patient's bedside—and easily obtain an adequate blood sample for single dose tests.



◇ Sterility cap

◇ Depth adjuster (preset at medium)

◇ Safer T-shape design



◇ Trigger



Color	Reorder#	Gauge	Penetration Depth
	05-2928A	28G	1.3/1.8/2.3mm
	05-2926A	26G	1.3/1.8/2.3mm
	05-2923A	23G	1.3/1.8/2.3mm
	05-2921A	21G	1.3/1.8/2.3mm

Newborn disease screening

Specialized Heel Incision Safety Lancets




BKAMED *Plus*

Heel Incision Safety Lancets

◇ One-handed blood drawing to increase efficiency

The one-piece safety lock can be released with one hand
Position the heel with one hand and push the button to fire



◇ Micro-blade with patented gentle launch mechanism and smooth blade path

Precise depth control and adequate blood volume without massaging reduces hemolysis, improves test accuracy and speeds the small incision healing




◇ Specialized product for newborns with assured quality and safety

Individually blister packaged and gamma radiation-sterilized
Blade retracts after use, preventing cross-contamination




Color	Reorder#	Incision Depth	Incision Length
	30-106514	0.65 mm	1.4mm
	30-101025	1.0 mm	2.5mm
	30-108575	0.85 mm	1.75mm
	30-102030	2.0 mm	3.0mm

Neoheel Heel Incision Safety Lancets



Shallow incision makes healing faster

Ordinary safety lancets



Deep incision makes healing slower

VS

Improving diabetes management
 Helping diabetics live a better life

Creating a "pain-free" blood sampling experience
 Enhancing blood sugar monitoring

The perfect partner for home-use blood sugar test strips



BKAMED

LDE4

Adjustable Lancing Device with Ejector

◇ Patented SSC (Soft Stop Control) technology absorbs shock and reduces pain



Soft-Stop Control

At the moment of puncture, the needle engages the SSC structure which absorbs the impact of the needle and stabilizes the puncture, through elastic deformation and energy conversion.



◇ **9 adjustable depth levels**

Easy to read settings



◇ Optional transparent AST cap



Patented inner cap cleaning technology

Traditional caps consist of inseparable inner and outer layers. During long-term lancet use, blood residue may collect on the inner cap surface, risking bacterial growth. The LDE4 features a cleaning setting (↕ Setting). At this setting, the outer cap can be removed, to enable easy cleaning of the inner cap.



Reorder#	Description	Color	Package
02-2201	9 Adjustable depths	Gray/White	Individually packaged 100 pcs/case



BKAMED Mino

Lancing Device with Ejector

- Innovative lancet ejector ring design**
Single pull removal, for efficient and convenient operation
- Compact and lightweight**
30% smaller in size
- Screw-on cap design**
No need to align, easy for the elderly to operate
- Easy to clean and disinfect**
The pen cap features an open cavity design, for easy cleaning and disinfection

BKAMED Disposable Blood Lancets

Soft & Soft Pro Softsure & Softsure Pro


- ◇ Comfort grip for needle placement
- ◇ International design benchmarked
- ◇ Well-established reputation for quality

Pro cap design
Safe disposal of fine needles

Softsure protection against fine needle deformation during loading

TPC pain-free technology absorbs shock and reduces vibration

Two models to choose from



The Mino lancing device's tip features a no-impact, no-shock design

The tip of a regular lancing device shakes upon impact



Mino has a lancet ejector

Mino SE does not have a lancet ejector

Product Model	Reorder#	Description	Color	Package
Mino	02-2601	6 Adjustable depths	Gray/White	Individually packaged 100 pcs/case
Mino SE	02-2701			



Color	Soft Reorder#	Color	Softsure Reorder#	Color	Soft Pro Reorder#	Color	Softsure Pro Reorder#	Gauge
	01-0132		01-1732		01-0232		01-1832	32G
	01-0130		01-1730		01-0230		01-1830	30G
	01-0128		01-1728		01-0228		01-1828	28G
	01-0126		01-1726		01-0226		01-1826	26G
	01-0121		01-1721		01-0221		01-1821	21G

